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T: +30 28210 63222 | WEB: www.mobiakcare.com

INTENDED USE:

Antidecubitus Cushions with air cells are designed to limit and delay the formation of decubitus. They are also intended for people with special needs or the elderly who own wheelchairs and have a sedentary lifestyle for a long time. They can also be used by users who do not have any symptoms and want to offer more comfort in their daily life. The Antidecubitus Cushions can be used by the user himself. When we place the cushion, it is important that the user gets up or is lifted up at regular intervals to limit the prolonged compression of the pelvis in the area of protruding bones.

FEATURES:

- Silky feel and softness for maximum soft tissue protection
- Cover with air permeability , easy to wash & dry
- Cover with fire retardant properties
- Ergonomic design, which reduces heat build-up due to the between the air cells gap
- Lightweight and easy to carry (with handle on the cover)
- Interconnected 3D air cells provide high relief & effective pressure distribution
- The internal airflow automatically adjusts while sitting, according to body structure and weight, for a better fit
- Easy adjustment of the pillow to the desired position (with special straps)
- Protective cushion coating of TPU (thermoplastic polyurethane) for easy cleaning
- Cushion codes **0804517**, **0804518**, **0804519** offer the opportunity of automatically adjustment of air flow, according to user needs based on innovative design.

BOX COMPONENTS

- One Cushion
- One Cover
- One Manual Air Pump, (**only for 0804514, 0804515, 0804516**)

SET-UP INSTRUCTIONS:

1. Place the cushion on the seat.
2. Open the valves and inflate the cushion with the manual pump (only for 0804514, 0804515, 0804516) or open the valve to inflate automatically – self inflate (only for 0804517, 0804518, 0804519)
3. Place the patient on the cushion
4. We recommend readjusting the cushion pressure every 10 days based on user preference.

CLEANING & WASHING INSTRUCTIONS

CUSHION

- Wash the cushion with antiseptic soap.
- Wipe the cushion using wet cloth: rinse in clear water then air dry it
- Do not use detergents or other solvents.
- Do not use wash machine or tumble dry.
- Do not use heat or steam autoclave.

COVER

- Wash at 30 °C, not ironing .
- For hygiene reasons, each cushion should be used by one patient only.

TECHNICAL SPECIFICATIONS:

CODE	DESCRIPTION	WITH	DIMENSIONS (WxLxH)	MAX LOAD
0804514	AIR CELL CUSHION "ANESIS I" 42x42x5cm TPU	AIR	42 x 42 x 5cm	160 Kg
0804515	AIR CELL CUSHION "ANESIS I" 40x46x5cm TPU	AIR	40 x 46 x 5cm	160 Kg
0804516	AIR CELL CUSHION "ANESIS I" 46x46x5cm TPU	AIR	46 x 46 x 5cm	160 Kg
0804517	AIR CELL SELF INFLATING CUSHION "ANESIS II" 42x42x5cm TPU	AIR	42 x 42 x 5cm	160 Kg
0804518	AIR CELL SELF INFLATING CUSHION "ANESIS II" 40x46x5cm TPU	AIR	40 x 46 x 5cm	160 Kg
0804519	AIR CELL SELF INFLATING CUSHION "ANESIS II" 46x46x5cm TPU	AIR	46 x 46 x 5cm	160 Kg

CODE	CUSHION COATING MATERIAL	COVER MATERIAL COMPOSITION	VALVES
0804514	TPU	UPPER PART: KNITTED STRETCH FABRIC WITH FIRE RETARDANT PROPERTIES	1
0804515	TPU		1
0804516	TPU		1
0804517	TPU	BOTTOM PART: PVC LEATHER WITH FIRE RETARDANT PROPERTIES	1
0804518	TPU		1
0804519	TPU		1

CAUTION(!)

- Do not smoke & do not use any flame near the cushion
- Do not expose the cushion to direct heat.
- If you need to store the cushion deflated first.
- Clean the product in weekly frequency.
- Protect the parts from scratches, cuts and punctures.
- In case of damage or malfunction, please contact your authorized and trained dealer.
- Do not use the product for a purpose not indicated in this manual
- MOBIAK S.A declines all responsibilities for any consequences resulting from an incorrect use of this product or use other than that mentioned in this manual.
- MOBIAK S.A reserves the right to make changes on the device and the information contained in this document without previous notice.

PROTECTING THE ENVIRONMENT

If your product needs replacement, or no longer serves you, consider the environment.

- 1) Do not dispose of your product along with the municipal waste (this is also the importance of the shown recycling symbol).
- 2) Contact your municipal authority to indicate the point of disposal of your product for recycling
- 3) You have to dispose of your useless product in the right place.

INCIDENT REPORTING:

Any serious incident that occurs in relation to use of the device must be reported to the manufacturer and a competent authority of the Member State in which the user and / or patient resides.

NORTON SCALE DETAILS

The Norton Scale was developed in the 1960s and is widely used to assess the risk for pressure ulcer in adult patients. The five subscale scores of the Norton Scale are added together for a total score that ranges from 5-20. A lower Norton score indicates higher levels of risk for pressure ulcer development.

SCORE	RISK FOR PRESSURE ULCER DEVELOPMENT
>18	LOW
18-14	MEDIUM
14-10	HIGH
< 10	VERY HIGH

Please be noted that the Norton Scale is a tool only to be used by trained healthcare professionals, and only as a supplement to a thorough patient state evaluation.

MBK CODE	SUITABLE FOR NORTON SCALE SCORE
0804514 0804515 0804516 0804517 0804518 0804519	19 – Low Risk

DECLARATION OF CONFORMITY:

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 745/2017 and where appropriate, the standards and legislation referred to.

WARRANTY:

The product has a warranty of good operation for two (2)* years from purchase date.

The warranty covers factory defects & does NOT cover defects resulting from misuse, poor maintenance, modification, overuse or non-compliance with the product instructions for use.

It also does NOT cover worn parts and whatever it is going to be damaged during use or over time.

Damages or defects caused by: natural disasters, unauthorized maintenance or repair, power supply problems (where provided), & transport companies are also not covered by the warranty.

*** six (6) months for the Manual Air Pump and Cover**

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Description of symbols that have been printed on the label, Master Carton and or User Manual

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